





**Catalent**<sup>®</sup>  
BIOLOGICS



## Best Practices, Strategies & Utilization of Novel Biological Responses for Robust Cell-Based Potency Assays


**PEDRO MORALES**  
DIRECTOR, BIOLOGICS ANALYTICAL

©2024 CATALENT, INC. ALL RIGHTS RESERVED.



more products. better treatments. reliably supplied.™

# INTEGRATED & STANDALONE CDMO SOLUTIONS TO ACCELERATE YOUR BIOLOGIC TO MARKET



# Biologics Analytical Services

Single Source for Integrated & Standalone Capabilities


## GMP Analytical Services

- Method development, transfer & optimization
- Phase-appropriate validation
- In-process, release & stability testing
- Binding & cell-based assays
- Extractables & leachables

## GLP Support for Clinical Studies

## Experience with broad classes of large molecules:

- Monoclonal, polyclonal & bispecific antibodies
- Bioconjugates & ADCs
- Oligonucleotides
- Recombinant proteins
- Fusion proteins
- Pegylated peptides
- Cell and gene therapies
- Aptamers
- Vaccines
- Oligosaccharides



**225+** scientists  
across the sites

**100,000+** ft<sup>3</sup> of  
stability chambers

**800+** assays/  
techniques offered


**300+** client programs  
supported

# AGENDA

- 1** Introduction & Considerations
- 2** Readouts
- 3** Strategy
- 4** Roadmap Visualization
- 5** Examples: Dose Response Curves
- 6** mRNA Transcription Assay Development & Qualification
- 7** Q&A



# Introducing Cell-Based Potency Assays



- **ICH Q6B** -> validated cell-based potency assays for commercial release & stability of large molecules for licensing in multiple countries.
- Must reflect the **MoA** of large molecules, be **stability indicating**, key for ensuring **Safety, Efficacy, Purity, Identity** and **Potency** of medicinal products.
- May be also utilized for **manufacture comparability assessments, extended characterization** and **formulation development**.
- ATMPs like CGT require a customized **“matrix”** panel of cell-based potency assays for manufacturing, purification and release decisions.

# General Aspects to Consider During Cell-Based Potency Method Development

## Cell Stimulation

- Direct analyte interaction
- Ligand interaction
- Serum/growth factor starvation
- Co-culture required
- ECM interaction necessary
- Others, such a biomechanical, hypoxia, etc.

## Biological Response

- The biological response should reflect the mechanism of action & be measurable
- Understand procedure to control for biological response drifts
- Start assay trending early to include relevant variables

## Curve Shapes

- Most biological responses can be analyzed with 4-PF curves, but straight-line analysis is also an option
- Ensure curve represents optimized method before validation

# AGENDA

- 1 Introduction & Considerations
- 2 **Readouts**
- 3 Strategy
- 4 Roadmap Visualization
- 5 Examples: Dose Response Curves
- 6 mRNA Transcription Assay Development & Qualification
- 7 Q&A



# Bioassays Readout - Where Biology Meets Chemistry!

Each of these assays use luminescent or fluorescent readouts to detect & quantify specific biochemical or cellular processes

Cyto-  
toxicity

Enzymatic

Reporter

Metabolic

Nucleic  
Acid

## Examples

- CytoTox-Glo™ assay measures a distinct protease activity associated with cytotoxicity
- alamarBlue™
- Caspase-Glo® assay measures caspase-3/7, enzymes that play a role in apoptosis
- cAMP-Glo™ assay measures cAMP, a second messenger important in many biological processes, by converting it to a luminescent signal
- CellTiter-Glo® assay measures ATP as an indicator of metabolically active cells
- alamarBlue™
- Quantitative PCR (qPCR) or droplet digital PCR (ddPCR), which are used to quantify DNA or RNA from cells



# AGENDA

- 1 Introduction & Considerations
- 2 Readouts
- 3 **Strategy**
- 4 Roadmap Visualization
- 5 Examples: Dose Response Curves
- 6 mRNA Transcription Assay Development & Qualification
- 7 Q&A




# Cell-Based Assay Strategy

## Part 1 – Cell culturing optimization

- Utilization of well characterized cells
- Evaluate continuous culture conversion to “thaw & use”
- Monitor critical parameters (reagents, cell viability, passage number, culture conditions & population doubling).
- GMP banks produced based on mock banks:
  - Master & working banks recommended
  - Generated & tested under protocol (QA & client approved)
  - Tested for mycoplasma, bioburden & functionality

## Part 2 – Assay

- Plate layout/uniformity (temperature/evaporation, etc.)
- Convert continuous culture to “thaw and use”
- Randomize distribution of replicates across plate
- Maximize throughput (horizontal vs vertical)
- Inclusion and trending assay control sample



## Pre-Validation Activities

- Performed under an agreed plan
- Involve multiple analysts/Plate readers

## Phase Specific Method Validation

- Under protocol with criteria per pre-validation
- Minimum of two analysts
- Based on client needs and minimum requirements of Catalent SOPs

# AGENDA


- 1 Introduction & Considerations
- 2 Readouts
- 3 Strategy
- 4 **Roadmap Visualization**
- 5 Examples: Dose Response Curves
- 6 mRNA Transcription Assay Development & Qualification
- 7 Q&A



# Roadmap to a Validated Cell-Based Assay

## Aspects to Consider

- Incubation times
- Critical reagents
- Training
- DS/DP scope
- Trending
- System suitability
- Locking steps
- Robustness
- Number of plates
- Critical steps
- Assay parameters
- True ACS lot
- Assay acceptance criteria
- Specifications
- Locking parameters
- Specificity




# AGENDA

- 1 Introduction & Considerations
- 2 Readouts
- 3 Strategy
- 4 Roadmap Visualization
- 5 **Examples: Dose Response Curves**
- 6 mRNA Transcription Assay Development & Qualification
- 7 Q&A



# Bioassay Dose Response Curves- Example of Optimal

- Drug dilution series generate biological responses
- Accuracy/Precision across range of method




Example of well-balanced curves, with defined upper and lower asymptotes, more than one point at the linear range, good parallelism and responses across the assay range

# Examples – Alternative Optimized Reagents

## Optimization of System Parameters

- Detection reagent choice: traditional versus optimized reagent



**Optimized Reagent**


**Traditional Reagent**

Trade-off between light intensity and duration of the signal using an optimized reagent

# Examples – Ligand Activity Optimization

## Optimization of System Parameters

- mAb binds to ligand and neutralizes activity (apoptosis in this example)
- Choice of ligand concentration critical for optimal dose response







# Examples – Incubation Conditions During Reactions

## Optimization of incubation with Luciferase substrate



- Identical dilutions used on two plates, one not cooled prior to substrate addition and the other cooled to RT for 15 min.
- Cooling may increase ambient O<sub>2</sub> solubility which is required during end point of reaction (Luciferase substrate requires ATP, oxygen and metallic cation).



# Examples – Incubation Time Effects

## Optimization of System Parameters


- Length of assay incubation dictated by drug MOA
- Comparison of curve shape across 3-day versus 4-day incubation



Extra day of incubation provides lesser steep curves, more points within hills-slopes, preventing “cobra shapes”, improved parallelism and goodness of fit

# Example - Cell Density Comparison vs Curve Shapes

○ RS □ ACS ▲ TS1 (200% RS) ◇ TS2 (50% RS)



Simple adjustment of cell density from 10,000 cells/well to 20,000 cells/well makes a difference in the response curves in terms of parallelisms and goodness of fit


# AGENDA

- 1 Introduction & Considerations
- 2 Readouts
- 3 Strategy
- 4 Roadmap Visualization
- 5 Examples: Dose Response Curves
- 6 **mRNA Transcription Assay  
Development & Qualification**
- 7 Q&A



# mRNA Transcription as Readout for Cell-base Potency


## Challenge



- The available prototype was a 6-well plate requiring RNA extraction prior cDNA synthesis
- Transcription regulation was required to recapitulate MOA of the product
- 6-well plate not feasible for GMP assay due to low throughput for all replicates, controls, and full curve
- However, a 96-well plate requiring RNA extraction will not provide enough RNA yield for downstream activities
- RNA extractions efficiencies accounts for variability, thus impacting the robustness of the assay


## Solution

- An optimized direct lysis method
  - Cell Lysis (~7 min total)
  - Reverse Transcription (~5 min)
  - Fast Real-Time PCR (~35 min)
- Faster direct cell lysis after stimulation of response
- No RNA extraction, no variability between extractions
- Higher throughput
- qPCR directly from lysates using TaqMan chemistries
- RNA can be easily converted to cDNA using the lysates and then amplified by PCR:



- TaqMan chemistries allows for multiplexing in same reaction (FAM vs VIC, for example), allowing well to well signal normalization



# Assay Schematical Overview



**Serum-Free supplemented?**

Yes

No



Dead

No

Yes

★ *Critical step identified – cells are thawed directly on serum-free medium to allow maximization of response to ligand.*

Cells previously treated with plasmid, ligand, or transgene are subjected to direct lysis




Cell lysate is subjected to reverse transcription to generate cDNA of transcript



The cDNA is then placed into a qPCR duplex qPCR reaction to assess cycle threshold ( $C_T$ ) of the target & reference transcripts

# General Assay Development

- Ligand and antagonist concentration titrated to find optimal induction concentration
- Fold-increase evaluated using  $\Delta\Delta C_T$  Analysis




Initial proof of concept with a limited number of data points.  
Ligand is tested at a fixed concentration and antagonist tested at titrated concentration to find optimal antagonistic working concentration


# mRNA Transcription - $\Delta C_T$ Analysis

## Principle

- Each well tested on qPCR is evaluated for cycle threshold ( $C_T$ ) for the target sequence & the reference sequence
- Calculate  $\Delta C_T$  (difference between  $C_T$  values of target & reference sequences) against concentration of the ligand or dilution factor of transgenic vector, and plot it
- Using  $\Delta C_T$  is advantageous because each well is normalized to each reaction. Small variations in total cDNA content among replicates becomes negligible



## Raw Data



- Calibrator gene = housekeeping gene with consistent transcriptional activity (no change with treatments)
- Amplicons (mRNA) that cross the cycle threshold earlier are more abundant or highly expressed after treatments
- Differential detection (FAM vs VIC) allows normalize well-to-well





# Qualification Plan

---

- Range of 5 concentrations/test samples tested (dilutional samples):
  - 50%, 70%, 100%, 140% and 200%
- On each plate, all samples (in 8-point serial dilution), run as pseudo-replicates:
  - Reference standard
  - 3 independent preparations of the same test sample
  - Assay control sample
- 3 plates run for each test sample (dilutional sample)
- 3 plates for 100% sample run by one analyst (for repeatability/intra-assay precision)
- 3 additional plates for 100% sample run by alternative analyst, and combined with first three 100% sample plates (for intermediate/inter-assay precision)

# Four Parameter Fit Analysis



4-P Fit:  $y = (A - D)/(1 + (x/C)^B) + D$ :


	A	B	C	D	Rel. Pot.
○ RefStd (RefStd: Concentration vs Values)	3.28	5.35	9.37	8.74	1
□ 50%TS (50PctTS: Concentration vs Values)	3.28	5.35	19.4	8.74	0.483
△ 100%TS (100PctTS: Concentration vs Values)	3.28	5.35	9.76	8.74	0.96
◇ 200%TS (200PctTS: Concentration vs Values)	3.28	5.35	5.43	8.74	1.73

Weighting: Fixed  
 PLA (Std. Curve: RefStd) Degrees of Freedom: parallel = 25 free = 16 non-parallel = 9  
 R<sup>2</sup> = 0.996 F-stat = 0.857 F-prob = 0.579

ΔCT is plotted against Sample Concentration.  
 Plots shown are 200% and 50% of Reference Standard and Control

# Qualification Results: Accuracy

## 50% Accuracy Level Test



	A	B	C	D	Rel. Pot.
○ RefStd (RefStd: Concentration vs Values)	3.22	4.46	7.15	8.56	1
■ 50%TS1 (50PctTS1: Concentration vs Values)	3.22	4.46	15.8	8.56	0.453
△ 50%TS2 (50PctTS2: Concentration vs Values)	3.22	4.46	16.1	8.56	0.445
◇ 50%TS3 (50PctTS3: Concentration vs Values)	3.22	4.46	15.3	8.56	0.467
● ACS (ACS: Concentration vs Values)	3.22	4.46	7.9	8.56	0.905

Weighting: Fixed

PLA (Std. Curve: RefStd) Degrees of Freedom: parallel = 32 free = 20 non-parallel = 12


R<sup>2</sup> = 0.997 F-stat = 0.826 F-prob = 0.624

Expected Potency (% Relative Potency)	Observed Potency (% Relative Potency)	Accuracy (% of expected)
50	46.7	93%
70	66.0	94%
100	100.2	100%
140	142.0	101%
200	175.1	88%

# Qualification Results: Precision

	Reportable Potency Values (Mean of 3 independent values from same plate)	Averaged Values across the three plates	Intra-Assay (Repeatability) Precision	Inter-Assay (Intermediate) Precision
<b>Analyst # 1</b>	104.3	100.2	3.3%	4.4%
	100.1			
	96.3			
<b>Analyst #2</b>	103.4	103.2	[Redacted]	
	109.3			
	96.9			

# Qualification Results: Specificity



4-P Fit:  $y = (A - D) / (1 + (x/C)^B) + D$ :

	<u>A</u>	<u>B</u>	<u>C</u>	<u>D</u>	<u>Rel. Pot.</u>
○ RS (RS: Concentration vs Values)	3.05	4.3	8.42	8.55	1
△ TS (TS - Rep 2: Concentration vs Values)	3.05	4.3	135	8.55	0.0624

# Comparison With Other Bioassay Readout Types

Assay Type	mRNA Transcriptional Assay by qPCR	Reporter Gene Assays
<b>Advantages</b>	Evaluates the mRNA expression of endogenous or exogenously introduced genes	Evaluates response to ligand directly (reporter)
	Because amplification it is more sensitive to differences in response	Options for off-the-shelf cell lines available
	Allows for analysis of transfected, transduced or edited gene activity	Evaluates functionally by product
<b>Limitations</b>	Analysis requires further transformation of raw data	Requires genetic engineering of cell line of interest

# Context For Cell and Gene Therapy Assays

Assay Type	ddPCR	mRNA Transcriptional Assay by qPCR	Flow Cytometry
<b>Advantages</b>	Assess absolute copy number transduced	Assess inserted or transferred gene construct functionality	Evaluates translational activity
	No standard curve required	Verification of cell transcriptional activity	Determines population of cells expressing transgene
<b>Limitations</b>	Low dynamic range	Requires normalization using a housekeeping gene	No amplification detectability limitations
	Unable to indicate if gene construct is inserted functionally	Multi-step assay	Difficult to standardize

# Conclusion



Assess response at mRNA transcription

No variability due to RNA extraction & purification

Well-by-well normalization results in better precision

Assay has been successfully qualified

Assay is accurate and precise and suitable for cGMP use

Additional tool for gene therapy with potential to conversion to ddPCR



# Context For Cell and Gene Therapy Assays

Fits into overall assay panel for products

Powerful when used in conjunction with orthogonal assays (e.g., flow cytometry, mass spectrometry)



Can evaluate performance of genes that impact metabolites (e.g., carbohydrates, proteoglycans, lipids)

Can be used when gene products are difficult to quantify and coupled with sequencing could be used to verify gene expression editing



# Catalent<sup>®</sup>

BIOLOGICS



## Q&A

©2024 CATALENT, INC. ALL RIGHTS RESERVED.



more products. better treatments. reliably supplied.™