


Modernizing the Platform Characterization Peptide Map for Accurate Assessment of Deamidation and Isomerization by LC-MS/MS


Lisa A. Marzilli, Elaine Stephens, Michelle English,
Shibu Philip* and Jason C Rouse

Analytical R&D, Pharmaceutical Sciences, Andover, MA, *Pearl River, NY



Antibody Considerations as a Drug Candidate

Typical parameters in selecting lead antibody:



Antibodies have a propensity to degrade:
manufacturing, storage and in vivo administration

Early Product Development Challenge

Determine
Degradation
Hotspots

Sequence
Liabilities or
CQAs?

Engineer to
Eliminate
Sequence
Liability?


- Oxidation: Met & Trp
 - **Deamidation: Asn and Gln**
 - **Isomerization: Asp**
 - Glycation: Lys
 - Truncations: Asp-Pro
 - Molecular modeling
 - Peptide mapping
- Does hotspot impact antigen binding?
 - Focus on complementarity determining regions (CDRs)
 - Constant region hotspots may be potential liabilities



Deamidation/Isomerization Scheme

- Asn deamidation reactions proceed more quickly at elevated pH and temperature
- N+1 residue (NG, NS, NT, and NN)

Geiger T. and Clarke S., J of Biol. Chem., (1987), 262 (2), 785.




- Asp isomerization reaction rate most impacted by lower pH and N+1 residue
- N+1 residue (DG, DS, DT, and DD)

Li et al, J. Pharm. Sci., (2013), 102 (3), 947.

The Ideal Peptide Map

- Robust and easy to perform
- Complete sequence coverage
- High quality separations (in a MS-friendly solvent)
- High quality UV chromatograms
- Minimal method-induced artifacts (short digestion times, low pH)
- Accurate targeted quantitation (Extracted Ion Chromatograms)



Developed in-house, automated software to highlight potential “hotspots”

Trastuzumab L Chain CDR Hotspots

1	DIQMTQSPSSLSASVGDRVTITCRASQDV <u>N</u>	30
31	<u>T</u> AVAWYQQKPGKAPKLLIYSASFLYSGVPS	60
61	RFSGSRSGTDFTLTISSLQPEDFATYYC <u>QQ</u>	90
91	<u>HYTTPPTFGQ</u> GTKVEIKR	

Trastuzumab H Chain CDR Hotspots

1	EVQLVESGGGLVQPGGSLRLSCAAS <u>GFNIK</u>	30
31	<u>DT</u> YIHWVRQAPGKGLEWVARIYPT <u>NG</u> YTRY	60
61	<u>ADSVK</u> GRFTISADTSKNTAYLQMNSLRAED	90
91	TAVYYCSR <u>WGGDGFYAM</u> DYWGQGTLLVTVSS	120
121		

- Automatically underlines CDR1, CDR2 and CDR3
- Highlights potential hotspots in blue



Case Study #1 Using Trastuzumab (Herceptin®) as a Benchmark for Peptide Map Method Development

Trastuzumab L Chain CDR Hotspots			Rel. Abundance
1	DIQMTQSPSSLSASVGDRVTITCRASQDV <u>N</u>	30	High level (>5%)
31	<u>T</u> AVAWYQQKPGKAPKLLIYSASFLYSGVPS	60	Low level (1-5%)
61	RFSGSRSGTDFTLTISSLQPEDFATYYC <u>QQ</u>	90	Potential (<1%)
91	<u>HYTTPPT</u> FGQGTKVEIKR		


Trastuzumab H Chain CDR Hotspots			Rel. Abundance
1	EVQLVESGGGLVQPGGSLRLSCAASGFNI <u>K</u>	30	High level (>5%)
31	<u>DT</u> YIHWVRQAPGKGLEWVARIYPT <u>NG</u> YTRY	60	Low level (1-5%)
61	<u>ADSVK</u> GRFTISADTSKNTAYLQMNSLRAED	90	Potential (<1%)
91	TAVYYCSR <u>W</u> GG <u>DG</u> FYAM <u>M</u> DYWGQGTLVTVSS	120	

Case Study #1 Using Trastuzumab (Herceptin®) as a Benchmark for Peptide Map Method Development

Trastuzumab L Chain CDR Hotspots			Rel. Abundance
1	DIQMTQSPSSLSASVGDRVTITCRASQDV <u>N</u>	30	High level (>5%)
31	<u>T</u> AVAWYQQKPGKAPKLLIYSASFLYSGVPS	60	Low level (1-5%)
61	RFSGSRSGTDFTLTISSLQPEDFATYYC <u>QQ</u>	90	Potential (<1%)
91	<u>HYTTPPT</u> FGQGTKVEIKR		


Trastuzumab H Chain CDR Hotspots			Rel. Abundance
1	EVQLVESGGGLVQPGGSLRLSCAAS <u>GF</u> NIK	30	High level (>5%)
31	<u>DTYIHWVRQ</u> APGKGLEWVARIYPT <u>NG</u> YTRY	60	Low level (1-5%)
61	<u>ADSVKGR</u> FTISADTSKNTAYLQMNSLRAED	90	Potential (<1%)
91	TAVYYCSRWGG <u>DG</u> FYAMDYWGQGTLLVTVSS	120	

Asn Deamidation and Asp Isomerization vs. Trypsin Digestion Time (Time Course) for Trastuzumab (Herceptin®)



Digestion Time Can Influence Level of Deamidation

Deamidation for Innovator Trastuzumab



- Time-course digestion facilitates the accurate measurement of deamidation at T=0
 - Digest for 30, 60 and 120 min
 - Report accurate deamidation levels


Literature	% Deamidation Level at N ⁵⁵ G (H-CDR2)
Genentech Harris 2001	~1 (CEX-HPLC)
Roche Sydow 2014	4 (pH 6.0 map)
Pfizer Results	
Innovator Trastuzumab	44 to 2

New Trypsin Peptide Map

Shorter Digestion Times and Time-Course

**Goal: minimize
method-related
artifacts**

**Peptide
Mapping**





- Time-course digestion facilitates the accurate measurement of deamidation at T=0
 - Digest for 30, 60 and 120 min
- Efficient reduction, alkylation and digestion with Lys-C/trypsin (Promega) (pH 8.2)
- Waters XSelect XP CSH C18 column(formic acid)
- Reproducible UV chromatographic profiles


- Analyzed on a Thermo Orbitrap Fusion™ Tribrid™ Mass Spectrometer
 - Unprecedented depth of analysis
 - High quality MS/MS data (CID, ETD)
 - Use Fusion for PTM discovery

Low-artifact Trypsin Digest Methods

Digest at pH 8.2



Digest at pH 6.0



Asn Deamidation and Asp Isomerization in Trastuzumab (Herceptin®) Using Various Methods



Site	Modification	Method	% Modification
N³⁰T (L-CDR1)	Deamidation	CEX-HPLC (UV) ¹	~15
		Trypsin pH 6.0 (overnight) map ²	11
		PFE-Trypsin pH 7.5 (overnight) map	12
		PFE-Lys-C/Trypsin pH 8.2 (T=0) map	10
N⁵⁵G (H-CDR2)	Deamidation	CEX-HPLC (UV) ¹	~1
		Trypsin pH 6.0 (overnight) map ²	4
		PFE-Trypsin pH 7.5 (overnight) map	44 
		PFE-Lys-C/Trypsin pH 8.2 (T=0) map	2
D¹⁰²G (H-CDR3)	Isomerization	CEX-HPLC (UV) ¹	~9
		Trypsin pH 6.0 (overnight) map ²	7
		PFE-Trypsin pH 7.5 (overnight) map	Not detected 
		PFE-Lys-C/Trypsin pH 8.2 (T=0) map	8

1. R.J. Harris, B. Kabakoff, F.D. et al, J of Chromat. B, Biomedical sciences and applications, 752, 233-245 (2001).


2. J.F. Sydow, F. Lipsmeier, V. et al, PLoS One, 9, e100736 (2014).

Efficient Separation and ETD Spectra Confirm ID of H12 (D¹⁰²G) and H12 (iso-D¹⁰²G)

H12 peptide (WGGD¹⁰²GFYAMDYWGQGTLVTVSSASTK)



Extracted Ion Chromatogram (EIC) for H12




• Waters XSelect XP CSH C18 column with formic acid can resolve isomerized peak




Isomerization of Asp-102 led to lower potency (Harris et al, 2001)

Diagnostic ion ($C_3 + 57$) Confirms Presence of isoAsp in Trastuzumab



Extracted Ion Chromatogram (EIC) for H12



- Waters XSelect XP CSH C18 column can resolve isomerized peak

Case Study #1

Peptide Map of Trastuzumab (Herceptin®)

- Trastuzumab (Herceptin®) is useful as a benchmark for peptide map method development
- Artifactual deamidation is a real possibility
 - Some motifs are sensitive to digestion time and pH (ex. N⁵⁵G in H-CDR2)
- Resolution of Asp isomerized peptides (D¹⁰²G vs isoD¹⁰²G) is challenging
 - Peptides have the same mass
 - Can only quantitate Asp/isoAsp if separated



Case Study #2: Deamidation in an IgG4 mAb

L Chain CDR Hotspots


1	XXXXXXXXXXXXXXXXXXXXXXXXXXXXRSSQSLA	30
31	NS Y NT FLSXXXXXXXXXXXXXXXXXXXXGISNRF	60
61	SXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX	90
91	XXXLQGT <u>HPY</u> TXXXXXXXXXX	

H Chain CDR Hotspots

1	XXXXXXXXXXXXXXXXXXXXXXXXXXXXGYRFT	30
31	NYW <u>I</u> HXXXXXXXXXXXXXXXXXXXXGINPG NN YATY	60
61	<u>RRKFQ</u> GXXXXXXXXXXXXXXXXXXXXXXXXXXXX	90
91	XXXXXXXXXEGYG NY <u>GAWFAY</u> XXXXXXXXXXXX	120
121		




Case Study #2: Deamidation in an IgG4 mAb



- Three distinct acidic isoforms (~15%)
- Acidic fractions collected and characterized for structure/function

Case Study #2: UV Chromatograms of Trypsin Peptide Maps




- CEX-HPLC Fraction A1 and A2 show an increase in H12 deamidation
 - Antigen binding impacted A1 charge variant


Case Study #2: New Peptide Map Column for Enhanced Chromatographic Separation

- IgG4 mAb: H chain deamidation in CDR3
- Localized deamidation site to: EGYGN¹⁰³YGAWFAYWGQGTLVTVSSASTK
- NY site not expected to deamidate (NG, NS, NT, NN)
- Causes loss of target binding

Previous Trypsin Map
Waters BEH C18, TFA



New Trypsin Map
Waters XSelect XP CSH C18, FA



Case Study #2: Deamidation of an IgG4 mAb

Site	Modification	Modification (total Asp and iso-Asp, %)	CEX Fraction
N³³S (L-CDR1)	Deamidation	3.2	A3
N⁵⁵N (H-CDR2)	Deamidation	0.1	- - -
N¹⁰³Y (H-CDR3)	Deamidation	7.0	A1 and A2

- New HPLC column allows for better separation and accurate quantitation of deamidation
- Detected unexpected N¹⁰³Y deamidation hot spot (beyond NG, NS, NT, and NN)
 - In process of updating software reflect new findings



Case Study #3

Forced Deamidation Study of IgG1 mAb

Stress Study

- Incubation at 25 °C and 37 °C in 100 mM sodium phosphate, pH 8.0
- T=0, 2 and 4 weeks

Peptide Map Method

- New Lys-C/trypsin map, pH 8.2, 30, 60 and 120 min
- Waters XSelect XP CSH C18 column, FA

MS Analysis

- LC-MS/MS and EICs
- Collision induced dissociation (CID)
- Targeted electron transfer dissociation (ETD)




Case Study #3: Forced Deamidation Study of IgG1 mAb

IgG1 L Chain CDR Hotspots			Rel. Abundance
1	DIVMTQTPLSLSVTPGQPASISCRSSQSLV	30	High level (>5%)
31	HS <u>NGNTFLYWYLQKPGQSPQLLIYRVS</u> NR <u>F</u>	60	Low level (1-5%)
61	<u>SGVPDR</u> FSGSGSGTDFTLKISRVEAEDVGV	90	Potential (<1%)
91	YYCF <u>QATHVPWTFGGG</u> TKVEIKR		

IgG1 H Chain CDR Hotspots			Rel. Abundance
1	QVQLVQSGAEVKKPGASVKVSCKASGYTFT	30	High level (>5%)
31	<u>GY</u> YIHWVRQAPGQGLEWMGWIYPG <u>NE</u> NTKY	60	Low level (1-5%)
61	<u>NERFKGRVTMTT</u> DTSTSTAYMELRSLRSDD	90	Potential (<1%)
91	TAVYYCARE <u>DG</u> SPYYAMDYWGQGTSVTVSS	120	


Deamidation of Light Chain (L2 peptide with CDR1)

25 SSQSLVHSN 33 GNTFLYWYLQKPGQSPQ 50 LLIYR 55





Deamidation of Heavy Chain (H5 with CDR2)

³⁹QAPGQ⁴³GLEWMGWIYPGN⁵⁵FNTK⁵⁹





CID Spectra Showing Deamidation at Asn-55

³⁹QAPGQGLEWMGWIYPGN⁵⁵FNTK⁵⁹



CID did not provide conclusive evidence for deamidation site

ETD Fragmentation Confirms Deamidation Site at Asn-55 ($N^{55}F$)



Asn Deamidation and Asp Isomerization Results for IgG1 mAb using New Low Artifact Trypsin Map

Site	Modification	Modification (total Asp and iso-Asp, %)		
		Unstressed	Stressed 4 weeks at 25 °C	Stressed 4 weeks at 37 °C
N³³G (L-CDR1)	Deamidation	2.8	12.8	43.0
N⁵⁵F (H-CDR2)	Deamidation	1.0	4.0	12.9
D¹⁰⁰G (H-CDR3)	Isomerization	0.7	1.0	1.7

- Important to analyze stressed material for hotspots (2.8 to 43.0%)
- Detected unexpected N⁵⁵F deamidation hot spot
- D¹⁰⁰G not a hotspot despite motif and location in H-CDR3

Summary

- **New trypsin peptide map-shorter digestion times, efficient enzyme, time-course elements and updated C18 column with formic acid**
 - Detects method-induced artifacts
 - Improved resolution of deamidated and isomerized peptides
 - Accurate quantitation of deamidation/isomerization via time course
 - Fragmentation (CID and ETD) determine sites of modification
 - ETD also distinguishes Asp vs iso-Asp
- **Deamidation and isomerization hotspots are sequence and location dependent and less predictable in proteins**
 - **Detected unexpected deamidation hot spots (beyond NG, NS, NT, and NN)**
 - N¹⁰³Y in heavy chain CDR3
 - N⁵⁵F in heavy chain CDR2
 - Low and constant level of Gln deamidation (Q⁵⁰L, Q⁴⁰G)

Summary of Hotspots

Case Study 1 Trastuzumab Originator	Case Study 2 IgG4 mAb Ref Material	Case Study 3 IgG1 mAb-stressed Ref Material - Stressed	Rel. Abundance
N³⁰T (L-CDR1)	N³³S (L-CDR1)	N³³G (L-CDR1) N³³G (L-CDR1)	High level (>5%)
N⁵⁵G (H-CDR2)	N⁵⁵N (H-CDR2)	N⁵⁵F (H-CDR2) N⁵⁵F (H-CDR2)	Low level (1-5%)
D¹⁰²G (H-CDR3)	N¹⁰³Y (H-CDR3)	D¹⁰⁰G (H-CDR3) D¹⁰⁰G (H-CDR3)	Potential (<1%)

Forced degradation studies help identify potential sequence liabilities



Acknowledgements

Jason Rouse-co-author
Elaine Stephens-co-author
Michelle English-co-author
Shibu Philip-co-author

Robert Kutlik (software)

Kristin Boggio
Heather DeGruttola
Jennifer Ide
Natalia Kozlova
Mellisa Ly
Andrew Saati
Roger Theberge
Matt Thompson
Ying Zhang

Tao He
Eric Sousa
Yan Zhang
Laura Lin

Bekah Ward
Neil Steinmeyer
Tom Porter
Meg Ruesch

